

Products and Order Codes

Standard E-Sacs

Product Code	D mm	L mm	V ml	Trocar mm	Number per box
EMP 70 INS	45	100	60	10	10
EMP 100 INS	55	160	225	10	10
EMP 130 INS	62	200	300	10;10/12	10
EMP 160 INS	102	250	850	10;10/12	5



Super E-Sacs

Product Code	D mm	L mm	V ml	Trocar mm	Number per box
EMP 100 MS	64	160	225	10/12	3
EMP 160 MN	102	230	900	12	3
EMP 180 MN	146	270	2000	15	3
EMP 200 MN	146	275	3000	15	3



Master E-Sacs

Product Code	D mm	L mm	V ml	Trocar mm	Number per box
EMP 80 WG	64	105	120	10	6
EMP 80 WG L	64	105	120	10	6
EMP 100 WG	64	150	225	10	6
EMP 100 WG L	64	150	225	10	6
EMP 120 WG	70	200	500	10	6
EMP 160 WG	108	240	900	15	3
EMP 180 WG	127	255	2000	15	3



Transfer Sleeves

Product Code	Trocar mm	Number per box
EMP 510 TS	5	6
EMP1215 TS	10/12	6

Sac Size	Surgical Procedure
70	Node biopsy, appendix, ectopic pregnancy, small gallbladder
80 80L	Standard gallbladder, small ovary or fibroid, adrenal, prostate
100 100L	Any gallbladder, ovary, small uterus, fibroid, prostate, adrenal, kidney
120 130	Ovary, fibroid, uterus, segmental bowel resection, larger adrenal, kidney
160	Ovarian cyst/tumour, large fibroid, uterus, kidney, spleen, bowel resection
180	Malignant kidney, large spleen, ovarian tumour, colon resection
200	Large tumours of ovary, kidney and spleen; hemi colectomy

